



Signet™ Y-SNP Identification System v2.0

Y-SNP Genotyping Core Kit

The Signet Y-SNP Identification System v2.0 comprises a core detection kit, a broad screening kit, and eleven individual genotyping kits. The Core Kit contains the Control DNAs, Exonuclease, 5x Labeling Mix, and Hybridization Solution. The screening and typing kits contain the specific 2X PCR Mixes and Bead Mixes necessary to complete the analysis.

All kits are for 100 Reactions.

Name	Cat. #	Name	Cat. #
Y-SNP Genotyping Core Kit	11710-100	Y-SNP KLMN Genotyping Kit	11717-100
Y-SNP Primary Screening Kit	11711-100	Y-SNP O1 Genotyping Kit	11718-100
Y-SNP AB Genotyping Kit	11712-100	Y-SNP O2 Genotyping Kit	11719-100
Y-SNP CD Genotyping Kit	11713-100	Y-SNP PQ Genotyping Kit	11720-100
Y-SNP E Genotyping Kit	11714-100	Y-SNP R1 Genotyping Kit	11721-100
Y-SNP FGHI Genotyping Kit	11715-100	Y-SNP R2 Genotyping Kit	11722-100
Y-SNP J Genotyping Kit	11716-100		

Note for customers familiar with Marligen Signet products: please pay attention to changes in the protocol!

An extensive Y-SNP bibliography and Y-SNP Haplotype Map are available from Marligen upon request from technical.service @Marligen.com

To Order: 866.464.4900 301.874.4990 customer.service@marligen.com

Important Information

The product you have received is authorized for laboratory research use only. The product has not been qualified or found safe and effective for any human or animal diagnostic or therapeutic application. Uses other than the labeled intended use may be a violation of applicable law.

Precautions

Warning: This product contains hazardous reagents. It is the end-user's responsibility to consult the applicable MSDS(s) before using this product. Disposal of waste organics, acids, bases, and radioactive materials must comply with all appropriate federal, state, and local regulations. If you have any questions concerning the hazards associated with this product, please call Marligen Biosciences, Inc at (301)874-4990.

Trademarks

Signet™ is a trademark of Marligen Biosciences, Inc. xMAP™ and Luminex® are trademarks of Luminex Corporation. MasterPlex™ is a trademark of MiraiBio, Inc. AmpliTaq Gold® is a trademark of Applied Biosystems Corporation. Platinum® is a trademark of Invitrogen Corporation.

Introduction to Signet™ Products

Signet products and reagents are based on xMAP™ Suspension Array Technology. DNA samples are amplified by PCR and labeled with a fluorescent tag. The fluorescently labeled products are then hybridized to an array of allele-specific oligonucleotides immobilized on latex beads. Each bead is identified by a unique spectral address (color). After hybridization, the Bead Mix is analyzed by the Luminex® 100 system that records the spectral address of the bead and the amount of fluorescent tag bound to it. Genotype reports can be generated by computer using a software package such as MasterPlex™ GT from MiraiBio.

Overview

The Marligen Y-SNP Identification System is designed to determine the genotype of human male Y chromosomes. Five simple steps comprise the assay:

1. Amplify the polymorphic Y-chromosome loci by multiplexed PCR.
2. Remove the PCR primers by incubating with Exonuclease I.
3. Label the PCR products by extension from a fluorescently-tagged oligonucleotide.
4. Hybridize the labeled DNA with allele-specific oligonucleotides in a Bead Mix array.
5. Read the hybridized microspheres in a Luminex 100 System

The Y-SNPs used to genotype the Y-chromosome are grouped into 12 multiplexes for highly efficient determination of a Y-chromosome genotype. The multiplexes are summarized in Table 1.

[2]

To Order: 866.464.4900 301.874.4990 customer.service@marligen.com



Table 1: Polymorphic loci in the Y-SNP Identification System and their multiplex assignments

Multiplex	A-R	AB	CD	E	FGHI	J
Multiplex description	Screening Multiplex	African, 168C	Asian, M168T/M89C/M96G	African, M Eastern, M96C	Eurasian M89T/ M9G/ M304A	M Eastern M304C
Variants	M9 M45 M89 M96 M122 M168 M175 M207 M304 M343 AmelXY	M31 M32 M42 M150 M146 M182 P4	M15 M38 M48 M55 M125 M130 M131 M145 M151 M210 M217	DYS391 M2 M33 M35 M58 M75 M78 M81 M123	M26 M52 M170 M201 M253 P15 P16 P37	M47 M67 M92 M172 M241 M267
Multiplex	KLMN	O1	O2	PQ	R1	R2
Multiplex description	Eurasian, M9C/ M45G/ M175+	SE-Asian, M175-/ M122T	SE-Asian, M175-/ M122C	Asian, Native American M45A/M207 A	Indo-European, M207G/M34 3C	Indo-European, M343A
Variants	M5 M11 M70 M178 M214 M231 SRY9138 Tat	M88 M95 M101 M103 M119 P31 SRY465	M7 M113 M117 M121 M134 M159 M164	M3 M19 M120 M143 M194 M199 M242 M323	M17 M56 M87 M124 M157 M173 SRY10831	M18 M37 M65 M126 M153 M269 P25 SRY-2627

The Multiplex A-R SNPs define major groupings of Y haplogroups. On samples of unknown origin, it is recommended to first test the samples with the screening multiplex A-R to determine the Y-haplogroups of the samples. The genotype can then be determined by subsequent testing with one of the Haplogroup-specific multiplexes.

Additional Materials Required

- Hot-start Taq DNA Polymerase (AmpliTaQ Gold® from Applied Biosystems or Platinum® Taq from Invitrogen)
- Taq Polymerase
- Water bath sonicator
- Vortex mixer
- Luminex 100 System with temperature controlled XY platform
- 96-well PCR plate and plate sealer or PCR tubes with caps
- Temperature-controlled block for 96-well plate set to 55°C (Optional)

[3]

To Order: 866.464.4900 301.874.4990 customer.service@marligen.com



Storage of Components

When you receive your Signet Y-SNP Identification System you may store the entire kit at -20°C. However, prior to first use of the kit the components must be separated and stored under different conditions.

Core Kit Components	Recommended Storage	Volume
9947A and 9948 DNA	-20°C	100 µl
Y-SNP Labeling Mix	-20°C, protect from light	500 µl
Exonuclease I	-20°C	100 µl
Hybridization Solution	Room Temperature	7 ml

Screening and Genotyping Kit Components	Recommended Storage	Volume
Y-SNP 2X PCR Mixes	-20°C, PCR set-up lab	1 ml
Y-SNP Bead Mixes	+4°C, protect from light	200 µl

General Precautions

Light Sensitivity

The Labeling Mix and Bead Mixes in this product contain fluorescent dyes that are light-sensitive. It is safe to use these items in normal laboratory fluorescent lighting. They should be covered and kept in a dark location when not in use. Aluminum foil is useful for creating temporary light shields. A drawer or darkened cabinet is also useful for short-term storage during use.

Bead Settling

Beads are 5.6 microns in diameter and are slightly more dense than water. They will settle gradually and should be mixed just prior to any pipetting step. Brief vortexing is sufficient for mixing suspensions that have settled for just a few minutes. Upon storage for several hours, small clumps of beads form that should be disrupted by sonication in a water bath. The Luminex 100 has a mixing function built into it that briefly injects a few microliters of sheath fluid prior to analyzing a sample. Therefore, disruption is not necessary during the plate reading.

Critical Parameters

- Use volumes precisely as indicated in the protocol
- Protect Labeling Mix and Bead Mixes from light
- Mix beads thoroughly to ensure a uniform suspension
- Mix PCR and Labeling Reagent by vortexing prior to use

[4]

To Order: 866.464.4900 301.874.4990 customer.service@marligen.com



PROTOCOL

I. Before beginning

The use of gloves and aerosol-resistant pipet tips is highly recommended. Keep pre-amplification and post-amplification reagents in separate rooms.

- 1) **Select the multiplexes to analyze.** On samples of unknown origin, it is recommended to first test the samples with the screening multiplex A-R to determine the Y-haplogroups of the samples. The genotype is then determined by subsequent testing with one of the Haplogroup-specific multiplexes. Different multiplexes may be run on the same PCR plate.
- 2) **Thaw the Y-SNP PCR Mixes.** It is very important to vortex the Y-SNP PCR mixes thoroughly for 5-10 seconds prior to each use. This dissolves any precipitate that may form during each freeze-thaw cycle.
- 3) **Check the Hybridization Solution for precipitate.** If a precipitate is present, heat the Hybridization at 55°C and mix until dissolved.
- 4) **Determine the number of reactions to be set up.** This should include a reagent blank (no DNA and no Taq Polymerase), a negative control (no DNA added but with Taq), and male DNA (9948) and female DNA (9947A) controls. Male DNA will yield a full profile, while female DNA will only give signal on the amelogenin X bead in Multiplex 1. Add 1 or 2 reactions to the total number of reactions to compensate for losses during pipetting. This will ensure that enough master mix is available for all samples in the experiment.

II. Amplify Y-chromosome DNA

- 1) Calculate the volumes of reagents needed for the master mix using the table below:

Component	Amount per Reaction	X Number Of Reactions	Final Volume
Y-SNP 2X PCR Mix	10 µl		
Taq DNA Polymerase (5 U/µl)	0.4 µl		
Control or test DNA	1-10 µl (1-10 ng)		
Deionized H2O	0-9 µl (To make 20 µl total volume)		

- 2) Prepare the master mixes by pipetting the volumes of everything except sample from the above table into 1.5 or 2 ml polypropylene tubes. Mix gently.

[5]

To Order: 866.464.4900 301.874.4990 customer.service@marligen.com



NOTE: The procedure requires a hot-start PCR protocol. We recommend use of Platinum Taq from Invitrogen.

- 3) Add 1-10 µl of each DNA sample (containing 1-8 ng) into a well or tube.
- 4) For the negative control, pipette nuclease-free water instead of DNA into the appropriate wells or tubes.
- 5) Pipette the appropriate volume of master mix into the wells or tubes.
- 6) Cap or seal the PCR plate, and place in a thermocycler.
- 7) Run the following Amplification Program

The duration of the initial 94°C step should be adjusted according to the instructions provided by the manufacturer of the Taq polymerase.

Initial Denaturation	94°C	Per manufacturer's instructions
Amplification: 30 cycles	94°C	30 sec
	58°C	20 sec
	72°C	90 sec
Final Extension and Hold	72°C	3 min
	5°C	Hold

Samples may be stored overnight at 4°C after the amplification step.

NOTE: If high sensitivity is not required, use ½ the recommended amount of Taq polymerase amount with the same reaction volume and/or scale the reaction to a final volume of 10 µl. Significant cost savings can be achieved with minimal reduction in sensitivity.

- 8) Remove primers with Exonuclease I

Add 1 µl Exonuclease I to each amplification reaction and mix thoroughly. Incubate at 37°C for 30 min. Inactivate the enzyme by heating at 80°C for 20 min. This step may be done conveniently in the thermal cycler.

III. Label the amplified DNA

A fluorescent tag is incorporated into the amplicons using Taq polymerase in a cycled primer extension reaction. Set these reactions up in a clean 96-well PCR plate. The total reaction volume of each well or tube is 25 µl.

NOTE: The Labeling Solution is light sensitive. Keep shielded from light.

- 1) Thaw the Y-SNP Labeling Mix. It is very important to vortex the Y-SNP Labeling Mix thoroughly for 5-10 seconds prior to each use to dissolve any precipitate that may form during each freeze-thaw cycle
- 2) Calculate the amount of Labeling Master Mix components that will be needed using the table below, adding one or two extra reactions to your calculation to ensure that

[6]

To Order: 866.464.4900 301.874.4990 customer.service@marligen.com



you have sufficient volume. The same Y-SNP Labeling Mix is used for all of the multiplexes.

Y-SNP Labeling Master Mix:

Component	Amount per Reaction	X Number of Reactions	Final Volume
diH ₂ O	18.75 µl		
Y-SNP Labeling Mix	5 µl		
Taq DNA Polymerase (5 U/µl)	0.25 µl		

- 3) Prepare the labeling master mix by pipetting the volumes of all components above into a sterile polypropylene tube. Mix gently.
- 4) Pipet 24 µl labeling master mix into the required number of sample and control wells or PCR tubes.
- 5) Add 1 µl PCR product to each well containing labeling master mix.
- 6) Cap or seal the PCR plate or tubes and place in a thermal cycler.
- 7) Run the following Labeling Program. Please note that the annealing temperature is 55°C rather than 58°C.

Initial Denaturation	94°C	Per manufacturer's instructions
Amplification: 40 cycles	94°C	30 sec
	55°C	30 sec
	72°C	30 sec
Final Extension and Hold	72°C	3 min
	5°C	Hold

Samples may be stored overnight at -20°C after labeling.

IV. Hybridize to Beads

The Signet™ Y-SNP Identification System contains Y-SNP Bead Mixes that match the PCR mixes. Make certain that a bead master mix is made for each multiplex or combination of multiplexes that were amplified. Different multiplexes may be analyzed in the same plate.

- 1) Warm up the Luminex 100 reader. Follow the manufacturer's instructions for preparing and calibrating the instrument. Using the Luminex Data Collector software, set the temperature of the XY platform to 55°C. Click on the Options icon followed by the Setup XY tab. Make certain that the "Heater Enabled" box is

checked, the Set Point is 55°C, and the aluminum 96-well plate holder is present in the XY platform.

NOTE: The Luminex 100 XY platform is required for this procedure. Make certain that the aluminum 96-well plate holder is installed and is pre-heated to 55°C.

NOTE: The Bead Mixes are light sensitive. Protect Bead Mixes from the light whenever they are not in use.

- 2) Add 35 µl Hybridization Solution to each well of the Costar 6509 96-well plate beginning at A1 and then continuing column-by-column until complete (B1-H1 followed by A2-H2, etc.).
- 3) Transfer the labeled PCR product into wells containing Hybridization Solution. Be certain to include the reagent blank controls and note their position. A reagent blank will be used as part of the machine set-up. Add the remaining samples column-by-column (A1-H1, followed by A2-H2, etc.) Place the plate into the pre-heated block of the Luminex 100 XY Platform and incubate at 55°C while preparing Hybridization Bead Mix.
- 4) Resuspend the Bead Mix by vortexing for 15-20 seconds. Sonicate the suspension in a water bath sonicator for 2 minutes to create a single-sphere suspension.
- 5) In a capped polypropylene tube prepare a Hybridization Bead Mix by adding 2 µl re-suspended Y-SNP Bead Mix to 35 µl of Hybridization Solution for each reaction. Make sufficient Hybridization Bead Mix for 1-2 more reactions than the number you will be setting up. Disperse the beads in the Hybridization Bead Mix by vortexing 15-20 seconds and then sonicating for 1 minute.
- 6) Add 37 µl of the Hybridization Bead Mix mixture to wells of the microtiter plate containing PCR product from the matching multiplex.
- 7) Hybridize at 55°C for 30 minutes. Protect samples from light.

V. Read in Luminex 100 System

Detailed instructions on hardware and software operations are provided in the Luminex 100 User's Manual. The following is an abbreviated protocol specific for the Signet Y-SNP Identification System. *For the bead map, please see the file available on Marligen's web site at www.marligen.com.*

- 1) The Luminex 100 should have been prepared before the hybridization step. A new session must be defined and the correct template loaded.
 - a) From the Luminex Data Collector screen, click Template on the toolbar.
 - b) A dialog box will appear. Select Signet Y-chromosome DNA II and click "OK".

[7]

To Order: 866.464.4900 301.874.4990 customer.service@marligen.com



[8]

To Order: 866.464.4900 301.874.4990 customer.service@marligen.com



NOTE: If Signet Y-chromosome DNA II does not appear on your list of templates, contact Marligen by e-mail at technical.service@marligen.com. The template and installation instructions will be sent to you via e-mail. Experienced users may use the bead map information from the Excel spreadsheet to key in their own template.

- c) A new dialog box will appear. Enter a name to identify this session on the New Folder line and click OK.
- d) The main session screen will appear. Click on the Setting button. A dialog box will appear. Enter the exact number of samples. The operator name may also be entered. The name of the session will appear in the Description box. All of the other boxes contain information specified in the template:

- e) Events: 150x (the number of beads) Total (i. e., 1500 total for J, 3300 total for A-R, etc.)
Sample size: 75 µl
Min events: 20
Flow rate: fast

Do not change any of these values. Click OK to save the completed settings.

- f) The main session screen should now contain a row for every sample. Descriptions of the samples may be entered by double-clicking on the individual cell and typing in the information.
- 2) Place samples in the XY tray. Fill the reservoir with sheath fluid and retract.
 - 3) Run one wash cycle.
 - 4) Check the Single box under the Start button. Make certain that the reagent blank is highlighted, and then click Start.
 - 5) The instrument will read the "Bead Only" sample in well A1.
 - 6) Use the data from the reagent blank to set the gate
 - a) Make sure that the X-axis of the histogram in the lower left corner is set to Doublet Discriminator
 - b) Right-click on the histogram. A pull-down menu appears
 - c) Scroll to Gate and follow the arrow to the right.
 - d) Click Create. Dotted, vertical red lines appear in the histogram frame.
 - e) With the mouse, drag each line to closely bracket the main peak. It is most important to keep the gate tight on the right-hand side, since this will eliminate microsphere doublets from the analysis.

[9]

To Order: 866.464.4900 301.874.4990 customer.service@marligen.com



- 7) Clear the Single box under the Start button. Make certain that the B1 row is highlighted, then click Start. The Luminex 100 will read the rest of the samples automatically.
- 8) After data collection is complete, select Save on the File menu. The data will be saved in the newly created session folder.

NOTE: If this step is omitted all data will be lost!

INTERPRETING RESULTS

The use of MasterPlex GT software from MiraiBio for data analysis is strongly recommended. MasterPlex GT conveniently automates the analysis of the output files produced by the Luminex Data Collector software. The main output file of Luminex Data Collect is called "Output.csv" and is located in the session folder that was named during session setup. Output.csv may also be imported into Microsoft Excel or similar spreadsheet programs for manual analysis.

The following parameters are recommended for data analysis with MasterPlex GT or for manual data analysis:

- 1) Use the PCR reaction(s) without Taq polymerase as the negative controls. This will ensure that the reagent background is subtracted from the data.
- 2) The minimum number of events (the threshold number of beads that must be recorded) should have a value of at least 20. This will avoid chatter generated by small numbers of stray microspheres.
- 3) Use the "relative intensity" for allele calling. To calculate the relative intensity of each allele, the MFI of the allele is divided by the sum of the MFI of both alleles:

$$RI_{A \text{ allele}} = MFI_{A \text{ allele}} / (MFI_{A \text{ allele}} + MFI_{B \text{ allele}})$$

- 4) Figure 1 shows the raw MFI data from 85 individuals at the M2 locus, collected on different days, machines, analysts, and DNA inputs. Although the MFI values vary considerably, the values cluster tightly around lines that correspond to a relative intensity of 88% for the G allele or a relative intensity of 68% for the A allele.

[10]

To Order: 866.464.4900 301.874.4990 customer.service@marligen.com






Figure 1. Raw MFI signals at locus M2. Although MFI levels can vary, relative intensity levels are nearly constant for each allele.

- 5) The intensity threshold should be set to at least 20. Beads with a median fluorescent intensity (MFI) less than 20 will not be reported as an allele. Samples with low intensity in all allelic forms (such as female DNA) will be reported as a null or no-call.
- 6) Figure 2 shows an example genotyping result for a test DNA that was tested with multiplex A-R. The results show that the sample is M9C/M45G/M175ins, therefore, this sample needs to be typed further with the KLMN multiplex.




Figure 2. Plot of Relative Intensity for the alleles of a Test DNA sample with the A-R multiplex. The sample type is M9C/M45G/M175ins, therefore, this sample needs to be typed further with the KLMN multiplex.

[11]

To Order: 866.464.4900 301.874.4990 customer.service@marligen.com



- 7) Figure 3 shows the genotyping result for the same test DNA sample tested with the KLMN multiplex. The results show that the sample is M214C/M231A/TatC/M178T, i. e., its haplotype is N3a.




Figure 3. Plot of Relative Intensity for the alleles of a Test DNA sample with the KLMN multiplex. The sample has been typed as M214C/M231A/TatC/M178T, and its haplotype is N3a.

Notes on some individual markers:

- Assays for some rare markers had been verified for the ancestral alleles only
- P25 is a multi-loci marker and will normally give a signal for both alleles. Users should call alleles by the C to A signal ratio, relying on the 9948 Control DNA, which is P25A (The haplotype of 9948 Control DNA is R1b3*). The signal ratio (C to A) for P25 marker is approx. 50:50 for A allele, or 60:40 for C allele.
- The P37 allele may generate a signal with female DNA samples due to cross-reactive sequences in the female genome. However, the results for P37 will be accurate when the sample comprises only male DNA.
- The assay for M17 will always develop strong signal for the *de*/allele. The signal for the *ins* allele is negligible for the *de*/type samples, or approx. 40% for the *ins* type samples.

Contact Us:

Marligen Biosciences, Inc.
2502 Urbana Pike
Ijamsville, MD 21754

301-874-4990
301-874-4993 (fax)
www.marligen.com

20109 rev. 081905

[12]

To Order: 866.464.4900 301.874.4990 customer.service@marligen.com

